



Xenon light source



Microfluidic
Reagent Rotor



Ready-to-Use reagents.
No warm-up period
required.

7 minutes

Sample-to-answer in
3 simple steps and
approximately 7-10 minutes

18 month


18-month
shelf life

27 tests

Up to 27 tests
in a single run


Clinical Chemistry⁺

Chemistry Analyzer Celercare® V5



MNCHIP

23 dedicated panel options covering 44 parameters



	# Comprehensive Diagnostic Profile (26+8) Comprehensive Diagnostic Profile (23+7) Comprehensive Diagnostic Profile (20+7) Comprehensive Diagnostic Profile (19+5) Health Check Plus Profile (15+5) Health Check Profile (14+5) Health Check Profile (13+4) Preanesthetic Plus Profile (9+2) Preanesthetic Profile (7+3) Critical Care Profile (13+4) Liver Profile (11+4) Liver & Kidney Profile (9+4) Kidney Profile (11+6) Electrolyte Profile (4+2) Diabetes Profile (4) Feline Inflammation Profile (10+3) Canine Inflammation Profile (7+2) Equine Profile (15+4) Avian & Reptile (12+4) Large Animal (10+3) Coagulation Profile (4) Triple tests Profile (3) Ammonia Profile																						
TP	TP	TP	TP	TP	TP	TP	TP	TP	TP		TP	TP	TP			TP	TP	TP	TP	TP		TP	
ALB	ALB	ALB	ALB	ALB	ALB	ALB	ALB		ALB		ALB	ALB	ALB			ALB	ALB	ALB	ALB	ALB		ALB	
GLO*	GLO*	GLO*	GLO*	GLO*	GLO*	GLO*	GLO*		GLO*		GLO*	GLO*	GLO*			GLO*	GLO*	GLO*	GLO*	GLO*		GLO*	
ALB/GLO*	A/G*	A/G*	A/G*	A/G*	A/G*	A/G*	A/G*		A/G*		A/G*	A/G*	A/G*			A/G*	A/G*	A/G*	A/G*	A/G*		A/G*	
ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT										ALT	
AST	AST	AST	AST	AST		AST		AST			AST	AST						AST	AST	AST			
AST/ALT*	AAR*	AAR*	AAR*	AAR*		AAR*		AAR*			AAR*	AAR*											
GGT	GGT	GGT	GGT	GGT							GGT	GGT				GGT		GGT		GGT			
ALP	ALP	ALP	ALP	ALP	ALP	ALP	ALP	ALP	ALP		ALP					ALP		ALP		ALP		ALP	
TBIL	TBIL	TBIL	TBIL	TBIL	TBIL	TBIL	TBIL				TBIL	TBIL						TBIL				TBIL	
DBIL	DBIL										DBIL												
IBIL*	IBIL*										IBIL*												
PT																					PT		
APTT																					APTT		
TT																					TT		
FIB																					FIB		
TBA	TBA	TBA	TBA	TBA	TBA						TBA					TBA			TBA				
CRE	CRE	CRE	CRE	CRE	CRE	CRE	CRE	CRE	CRE	CRE		CRE	CRE			CRE		CRE				CRE	
BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN			BUN		BUN		BUN		BUN	
BUN/CRE*	B/C*	B/C*	B/C*	B/C*	B/C*	B/C*	B/C*	B/C*	B/C*	B/C*		B/C*	B/C*			B/C*		B/C*				B/C*	
UA	UA																	UA					
CK	CK	CK		CK		CK	CK	CK										CK	CK	CK		CK	
LDH								LDH										LDH					
AMY	AMY	AMY	AMY	AMY	AMY	AMY	AMY									AMY	AMY					AMY	
LPS	LPS	LPS	LPS													LPS	LPS						
GLU	GLU	GLU	GLU	GLU	GLU	GLU	GLU	GLU	GLU	GLU		GLU	GLU		GLU		GLU	GLU	GLU			GLU	
FRU	FRU														FRU								
TG				TG											TG								
CHOL	CHOL	CHOL	CHOL	CHOL	CHOL	CHOL	CHOL				CHOL				CHOL							CHOL	
K+	K+	K+	K+		K+					K+			K+	K+				K+	K+				
Na+	Na+	Na+	Na+		Na+					Na+		Na+	Na+					Na+	Na+				
Na+/K+	Na+/K+	Na+/K+	Na+/K+		Na+/K+					Na+/K+		Na+/K+	Na+/K+					Na+/K+	Na+/K+				
Cl-	Cl-	Cl-	Cl-							Cl-		Cl-	Cl-						Cl-				
Ca	Ca	Ca	Ca	Ca	Ca	Ca	Ca			Ca		Ca				Ca	Ca	Ca	Ca	Ca		Ca	
P	P	P	P	P	P	P	P			P		P						P	P			P	
CaxP*	CaxP*	CaxP*	CaxP*	CaxP*	CaxP*	CaxP*	CaxP*			CaxP*		CaxP*						CaxP*	CaxP*			CaxP*	
Mg	Mg	Mg		Mg						Mg										Mg			
tCO ₂	tCO ₂	tCO ₂	tCO ₂	tCO ₂						tCO ₂		tCO ₂	tCO ₂					tCO ₂					
NH ₃																						NH ₃	
LAC	LAC	LAC								LAC													
pH										pH													
AG*	AG*	AG*	AG*							AG*		AG*	AG*										
fSAA																fSAA							
cCRP																	cCRP						

* Calculated test value # Coming soon