



Xenon light source



Microfluidic
Reagent Rotor



Ready-to-Use reagents.
No warm-up period
required.

7 minutes

Sample-to-answer in
3 simple steps and
approximately 7-10 minutes

18 month

18-month
shelf life

27 tests

Up to 27 tests
in a single run

Clinical Chemistry⁺


Chemistry Analyzer Pointcare[®] V3



MNCHIP

Pointcare® V3 Reagent Discs

23 dedicated panel options covering 44 parameters



Comprehensive Diagnostic Profile (26+8)
 Comprehensive Diagnostic Profile (23+7)
 Comprehensive Diagnostic Profile (20+7)
 Comprehensive Diagnostic Profile (19+5)
 Health Check Plus Profile (15+5)
 Health Check Profile (14+5)
 Health Check Profile (13+4)
 Preanesthetic Plus Profile (9+2)
 Preanesthetic Profile (7+3)
 Critical Care Profile (13+4)
 Liver Profile (11+4)
 Liver & Kidney Profile (9+4)
 Kidney Profile (11+6)
 Electrolyte Profile (4+2)
 Diabetes Profile (4)
 Feline Inflammation Profile (10+3)
 Canine Inflammation Profile (7+2)
 Equine Profile (15+4)
 Avian & Reptile (12+4)
 Large Animal (10+3)
 Coagulation Profile (4)
 Triple tests Profile (3)
 Ammonia Profile

TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP	TP
ALB	ALB	ALB	ALB	ALB	ALB	ALB	ALB		ALB		ALB	ALB	ALB				ALB	ALB	ALB	ALB	ALB		ALB
GLO*	GLO*	GLO*	GLO*	GLO*	GLO*	GLO*	GLO*		GLO*		GLO*	GLO*	GLO*				GLO*	GLO*	GLO*	GLO*	GLO*		GLO*
ALB/GLO*	A/G*	A/G*	A/G*	A/G*	A/G*	A/G*	A/G*		A/G*		A/G*	A/G*	A/G*				A/G*	A/G*	A/G*	A/G*	A/G*		A/G*
ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT	ALT											ALT
AST	AST	AST	AST	AST		AST		AST			AST	AST							AST	AST	AST		
AST/ALT*	AAR*	AAR*	AAR*	AAR*		AAR*		AAR*			AAR*	AAR*											
GGT	GGT	GGT	GGT	GGT							GGT	GGT					GGT		GGT		GGT		
ALP	ALP	ALP	ALP	ALP	ALP	ALP	ALP	ALP	ALP		ALP						ALP		ALP		ALP		ALP
TBIL	TBIL	TBIL	TBIL	TBIL	TBIL	TBIL	TBIL				TBIL	TBIL							TBIL				TBIL
DBIL	DBIL										DBIL												
IBIL*	IBIL*										IBIL*												
PT																						PT	
APTT																						APTT	
TT																						TT	
FIB																						FIB	
TBA	TBA	TBA	TBA	TBA	TBA						TBA						TBA			TBA			
CRE	CRE	CRE	CRE	CRE	CRE	CRE	CRE	CRE	CRE	CRE		CRE	CRE				CRE		CRE				CRE
BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN	BUN				BUN		BUN		BUN		BUN
BUN/CRE*	B/C*	B/C*	B/C*	B/C*	B/C*	B/C*	B/C*	B/C*	B/C*	B/C*		B/C*	B/C*				B/C*		B/C*				B/C*
UA	UA																			UA			
CK	CK	CK		CK		CK	CK	CK											CK	CK	CK		CK
LDH								LDH											LDH				
AMY	AMY	AMY	AMY	AMY	AMY	AMY	AMY										AMY	AMY					AMY
LPS	LPS	LPS	LPS														LPS	LPS					
GLU	GLU	GLU	GLU	GLU	GLU	GLU	GLU	GLU	GLU	GLU		GLU	GLU			GLU		GLU	GLU	GLU			GLU
FRU	FRU															FRU							
TG				TG												TG							
CHOL	CHOL	CHOL	CHOL	CHOL	CHOL	CHOL	CHOL				CHOL					CHOL							CHOL
K ⁺	K ⁺	K ⁺	K ⁺		K ⁺					K ⁺			K ⁺	K ⁺					K ⁺	K ⁺			
Na ⁺	Na ⁺	Na ⁺	Na ⁺		Na ⁺					Na ⁺			Na ⁺	Na ⁺					Na ⁺	Na ⁺			
Na ⁺ /K ⁺	Na ⁺ /K ⁺	Na ⁺ /K ⁺	Na ⁺ /K ⁺		Na ⁺ /K ⁺					Na ⁺ /K ⁺			Na ⁺ /K ⁺	Na ⁺ /K ⁺					Na ⁺ /K ⁺	Na ⁺ /K ⁺			
Cl ⁻	Cl ⁻	Cl ⁻	Cl ⁻							Cl ⁻			Cl ⁻	Cl ⁻						Cl ⁻			
Ca	Ca	Ca	Ca	Ca	Ca	Ca	Ca			Ca			Ca				Ca	Ca	Ca	Ca	Ca		Ca
P	P	P	P	P	P	P	P			P			P							P	P		P
CaxP*	CaxP*	CaxP*	CaxP*	CaxP*	CaxP*	CaxP*	CaxP*			CaxP*			CaxP*						CaxP*	CaxP*			CaxP*
Mg	Mg	Mg		Mg						Mg											Mg		
tCO ₂	tCO ₂	tCO ₂	tCO ₂	tCO ₂						tCO ₂			tCO ₂	tCO ₂					tCO ₂				
NH ₃																							NH ₃
LAC	LAC	LAC								LAC													
pH										pH													
AG*	AG*	AG*	AG*							AG*			AG*	AG*									
fSAA																	fSAA						
cCRP																		cCRP					

* Calculated test value # Coming soon